

# Improving Sperm Resistance

Viability following incubation at 37°C in function of time

TIME	CONTROL	maXipig®	SIGNIFICANCE
0 min.	93,5 ±2,0%	<b>92,1 ±1,6%</b>	NS
15 min.	92,8 ±2,4%	<b>93,0 ±1,1%</b>	NS
30 min	92,1 ±1,9%	<b>93,0 ±2,1%</b>	NS
60 min	66,8 ±1,0%	<b>90,5 ±1,9%</b>	*
90 min	60,3 ±1,2%	<b>94,6 ±2,3%</b>	*

Source: Yeste M. et al. **Specific LED-based red light photo-stimulation procedures improve overall sperm function and reproductive performance of boar ejaculates.** Scientific Reports 2016 Mar 2;6:22569

The logo for maXipig, featuring the word "maXipig" in a teal, sans-serif font. The letter "X" is significantly larger and bolder than the other letters. A stylized, teal-colored "6" is positioned to the right of the word, partially overlapping the "g".

# Improving Sperm Resistance

Acrosome integrity following incubation at 37°C in function of time.

TIME	CONTROL	maXipig®	SIGNIFICANCE
0 min.	95,0 ±2,1%	<b>95,8 ±2,3%</b>	NS
15 min.	94,1 ±2,0%	<b>95,8 ±2,3%</b>	NS
30 min	93,7 ±1,5%	<b>95,0 ±1,9%</b>	NS
60 min	85,1 ±1,1%	<b>96,3 ±2,4%</b>	*
90 min	70,7 ±1,1%	<b>92,8 ±2,5%</b>	*

Source: Yeste M. et al. **Specific LED-based red light photo-stimulation procedures improve overall sperm function and reproductive performance of boar ejaculates.** Scientific Reports 2016 Mar 2;6:22569



# Improving Sperm Resistance

Total Motility following incubation at 37°C in function of time.

Time	CONTROL	maXipig®	SIGNIFICANCE
0 min.	95,2 ±2,2%	<b>95,4 ±2,1%</b>	NS
15 min.	62,1 ±6,4%	<b>67,9 ±7,8%</b>	NS
30 min	62,8 ±5,1%	<b>63,0 ±6,1%</b>	NS
60 min	45,9 ±6,2%	<b>67,9 ±6,0%</b>	*
90 min	41,8 ±6,7%	<b>79,6 ±7,0%</b>	*

Source: Yeste M. et al. **Specific LED-based red light photo-stimulation procedures improve overall sperm function and reproductive performance of boar ejaculates.** Scientific Reports 2016 Mar 2;6:22569



# Fertility and Prolificacy

“In vivo” Fertility parameters of boar sperm ejaculates subjected to a previous photo-stimulation procedure

Treatment	N	Farrowing rate (%)	Total piglets at parturition	Live-born piglets at parturition
Control	800	83,7	13,5±0,2	12,7±0,2
maXipig®	520	88,1*	14,9±0,3*	14,0±0,2*

N = Number of sows

Source: Yeste M. et al. **Specific LED-based red light photo-stimulation procedures improve overall sperm function and reproductive performance of boar ejaculates.** Scientific Reports 2016 Mar 2;6:22569



# Fertility

“In vivo” Farrowing rate of boar sperm ejaculates subjected to a previous photo-stimulation procedure. Extended tests\*.

Treatment	N	Farrowing rate (%)
Control	6833	90,0
Photo-stimulated	3044	92,3

N = Number of sows

\* 2 years tests carried in 9 farms. Data from Ge Pork – alea!

maXipig<sup>6</sup>




alea!

# Fertility

Range	Nº sows CONTROL	Nº sows LED	% Fertility CONTROL	% Fertility LED	LED batch improvement
≤80% fertility	118	116	77,119	85,345	+8,226
80-85% fertility	1338	191	83,931	88,482	+4,551
85-90% fertility	1865	1032	89,383	91,472	+2,089
90-95% fertility	1958	1407	90,500	93,105	+2,605
≥95% fertility	1554	298	96,396	96,979	+0,583


difference between range groups of fertility index farms



\* 2 years tests carried in 9 farms. Data from Ge Pork – alea!

# Fertility

Fertility	n° sows CONTROL	n° CONTROL repeated sows	CONTROL Fertility	n° sows LED	n° LED repeated sows	LED Fertility	LED batch improvement
Summer 2014	328	53	83,841	207	23	88,889	+5,048
Autumn 2014	411	42	89,781	278	22	92,086	+2,305
Winter 14-15	372	26	93,011	213	11	94,836	+1,825
Spring 2015	483	39	91,925	323	15	95,356	+3,431
Summer 2015	2141	187	91,266	910	73	91,978	+0,712
Autumn 2015	2836	341	87,976	856	66	92,290	+4,314



\* 2 years tests carried in 9 farms. Data from Ge Pork – alea!

# Prolificacy

“In vivo” Fertility parameters of boar sperm ejaculates subjected to a previous photo-stimulation procedure

Treatment	N	Total piglets at parturition	Live-born piglets at parturition
Control	3822	13,22	12,15
<b>maXipig®</b>	<b>2179</b>	<b>14,00</b>	<b>12,85</b>

N = Number of births

\* 2 years tests carried in 9 farms. Data from Ge Pork – alea!

The logo for maXipig, featuring the word "maXipig" in a teal, sans-serif font, with a stylized piglet icon integrated into the letter "g".



# Farms

“*in vivo*” test in productive farms

Farm	Start date	Finish date	N° sows	Genetic line	Semen supplier	AI type	Boar breed
5 *	02/08/2014	31/12/2015	1320	A	a	60ml/2 billion sp	Duroc / Pietrain
2	28/04/2015	31/12/2015	750	A	a	60ml/2 billion sp	Landrace / Largewhite / Duroc / Pietrain
3	01/07/2015	31/12/2015	2200	A	b	90ml/3 billion sp	Pietrain
4	02/07/2015	31/12/2015	1500	B	c	45ml/1,5 billion sp	Pietrain
7	02/07/2015	31/12/2015	1500	C/D	a	45ml/1,5 billion sp	Duroc
6	06/07/2015	31/12/2015	500	A	a	45ml/1,5 billion sp	Pietrain
9	10/07/2015	31/12/2015	750	B	a	45ml/1,5 billion sp	Pietrain
8	16/07/2015	31/12/2015	1000	A	a	90ml/3 billion sp	Pietrain
10	02/09/2015	31/12/2015	800	A	a	45ml/1,5 billion sp	Pietrain

\* First “*in vivo*” test farm

maXipig<sup>6</sup>

