iPeak® reader 4,3"

Instructions for Use

Photometric immunoassay reader

Software version 4.3 and above

 EN / ES

Doc. 50008314-14

Revision date: May 2026

All Rights Reserved © 2026

Contents

Contents 3

ENGLISH 4

ESPAÑOL 69

Technical Support Contact Information 137

ENGLISH

Indication of the Introduced Modifications.....	6
Glossary	7
1.Introduction	8
1.1. Intended Purpose / Use of the iPeak® reader 4,3"	8
1.2. Intended Users & Intended Environment For Use	8
1.3. User Management	8
1.4. Device description	8
1.5. Contraindications	9
1.6. General information	9
2.Safety Information	12
2.1. Proper Use	12
2.2. Electrical Safety	13
2.3. Environment	15
2.4. Biological Safety	15
2.5. Chemicals	16
2.6. Maintenance Safety	17
2.7. Waste Disposal	17
2.8. Glossary of symbols	17
2.9. Cybersecurity	19
3.General Description	20
3.1. Device Overview	20
3.2. Accessories of the device	20
4.Installation	22
4.1. Unpack the device	22
4.2. Site requirements	22
4.3. Power cable connection	22
4.4. Cassette Requirements	22
5.Operating Procedures.....	23
5.1. Login screen	23
5.2. Home screen overview	24
5.3. Cassette Processing	24
5.4. System Menu	31
5.5. Side Menu	41
5.6. Tools Menu	49
5.7. Software update	52
5.8. Browser Generic Operations	52
6.Troubleshooting.....	54

6.1. Error and troubleshooting list	54
6.2. Remote Support	56
7.Maintenance.....	57
7.1. Cleaning Procedure	57
7.2. Service Life and Preventive Maintenance	57
APPENDIX A Ordering codes.....	58
APPENDIX B Technical data	60
APPENDIX C Warranty	62
APPENDIX D Waste Electrical and Electronic Equipment (WEEE)	63
APPENDIX E EU Declaration of Conformity	64
APPENDIX F FCC Statement.....	65
APPENDIX G UK Declaration of Conformity.....	66
APPENDIX H RoHS Statement	67
APPENDIX I Applied Standards	68

Indication of the Introduced Modifications

If it's the first time that you use the iPeak® reader 4,3" , please read the entire document carefully.

If it's not the first time that you use iPeak® reader 4,3" , you can read the indication of the introduced modifications in the table below.

Revision	Description of change
14	Update section 1.4.1, 1.6.1, 2.2, 3.2.4, 5.5.2, Appendix A, B, G, I

Glossary

LFR	Lateral Flow Reader, it is a device such as iPeak® reader 4,3" that reads and processes immunochromatographic lateral flow tests.
Administrator	This role holds the highest operational authority and has exclusive access to restricted functions of the LFR.
Manager	This role holds the second highest operational authority can manage users, printers, and network setups of the LFR.
Operator	This role holds the lowest operational authority. It is a person designated by the health organization to use the LFR for the intended use.
LFR Configuration	Set of options that determine the behavior and functionalities of the LFR. Most of them are only available to the Administrator user.
IUL's Partner	It is the institution, company, or other legal entity responsible to provide the LFR together with the cassettes to the health organization. Moreover, IUL's partner defines the necessary LFR's configuration for the intended use environment using the LFR's Administrator user mode to update the LFR's configuration when it is necessary.
Health Organization	It is an institution, company, or other legal entity responsible to operate the LFR together with the cassettes to be used in medical analysis of different nature.
Cassette	It is a device that consists of lateral flow tests and a housing that contains them.
Cassette Processing	The process that analyses lateral flow tests, from inserting a cassette into an LFR, until the results are shown on the screen. This process includes the scan of the cassette.
Cassette Model	It is the cassette housing, with a specific geometry that is adjusted in the slot of a tailored drawer.
Cassette Type	It is the set of a cassette model and specific lateral flow tests contained in its strips.
Cassette Settings	It is a set of parameters that enable the LFR to process and measure a specific batch of a specific cassette type. Also referred to as "Test Name", "Cassette Configuration" and "Cassette Config."
Data Matrix	Barcode that contains the full information of one cassette settings. Also referred to as "Configuration Barcode".
Small-Data Matrix	Barcode that contains the cassette settings code and the batch ID. It is used to confirm that the current cassette matches with current cassette settings. Also referred to as "Small Barcode" and "Confirm Cassette".

1. INTRODUCTION

Thank you for choosing the iPeak® reader 4,3" . We are confident that this instrument will become an integral part of your laboratory.

Before using the iPeak® reader 4,3" , it is essential that you read these Instructions for Use carefully. Following the instructions and safety information in these Instructions for Use will ensure safe operation and maintain the system in a safe condition.

In case a serious incident occurred in relation to the device, it shall be reported to IUL and the competent authority of the Member State in which the user and/or the patient is established.

1.1. Intended Purpose / Use of the iPeak® reader 4,3"

The iPeak® reader 4,3" is an instrument specifically to be used to provide quantitative, semi-quantitative, and/or qualitative in-vitro determination of photometric immunochromatographic test defined and marketed by IUL's partner.

The iPeak® reader 4,3" is intended to be used only in combination with lateral flow (LF) tests indicated for use with the iPeak® reader 4,3" , and only for applications that are described in the respective handbooks of our reader associates.

1.2. Intended Users & Intended Environment For Use

The iPeak® reader 4,3" is intended for professional use in clinical laboratories, hospitals, and outpatient settings (e.g., physician's offices). It is not intended for use in intensive care units or surgical environments, unless all specific hygiene and patient safety requirements of such settings are met.

1.3. User Management

The iPeak® reader 4,3" can be used with three different roles:

- **Operator:** have access to the features needed when operating.
- **Manager:** have access to features needed when operating, plus the network and connectivity features and can also delete scanned results.
- **Administrator:** have unlimited access to every feature of the iPeak® reader 4,3" .

1.4. Device description

1.4.1. Lateral flow reader

The iPeak® reader 4,3" is a portable, point-of-care (PoC) diagnostic instrument designed for the objective interpretation of colorimetric immunochromatographic lateral flow assays. The device has a compact footprint suitable for use in laboratory and near-patient testing environments.

The instrument provides qualitative, semi-quantitative, and/or quantitative analysis of compatible lateral flow test strips and cassettes. Results are generated through optical detection and analysis of the test and control lines present on the assay device.

The reader is designed to accommodate strips and cassettes of different sizes and formats as specified by the test manufacturer.

The instrument incorporates integrated data management functionality. Test results may be stored within the device and exported through multiple data transfer options. The reader supports connectivity to a Laboratory Information Management System (LIMS) server for electronic transmission of results.

The system includes user management functionality with configurable access levels to control system access. An electronic audit trail records user activities and system events, ensuring traceability of actions performed on the device, including test execution, result modification (if applicable), and configuration changes.

1.4.2. Development Configuration

Test-specific analysis settings are defined and validated using a dedicated configuration tool available to authorized partners. Once configured, these settings are uploaded to the Lateral Flow Reader and applied during test analysis.

1.5. Contraindications

Not applicable.

1.6. General information

1.6.1. Scope of delivery

The iPeak® reader 4,3" main unit is supplied together with the essential components required for its intended operation.

The adjustment drawer (Cat. 900013732) is provided as an accessory intended exclusively for corrective camera readjustment procedures in case of detected deviation of capture parameters. The adjustment drawer is not required for routine operation and does not form part of the essential intended purpose of the main unit. See this section on Adjustment drawer for more information.

1. iPeak® reader 4,3" main unit (Cat. 100033000).
2. Drawer (Customized for each IUL's Partner Cat. 9000xxxxx).
3. Adjustment drawer (Cat. 900013732)
4. External power supply unit (Cat. 900020088 or 900020089).
5. Power cord (Cat. 900020077, Cat. 900020079, Cat. 900020080, Cat. 900020090, or Cat. 90002264).
6. Lateral Flow Reader Datamatrix (Cat. 50009100)

See 'Ordering Codes' Appendix for more information.

1.6.2. Technical Assistance

At IUL we pride ourselves on the quality and availability of our technical support. Our Technical Service Department is staffed by experienced technicians with extensive practical and theoretical expertise in the use of IUL products. If you have any questions or experience any difficulties regarding the iPeak® reader 4,3" or IUL products in general, do not hesitate to contact us.

IUL customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other customers as well as to the researchers at IUL. We, therefore, encourage you to contact us if you have any suggestions about product performance or new application and techniques.

For technical assistance, contact the Technical Service Department.

1.6.3. Policy Statement

It is the policy of IUL to improve products as new techniques and components become available. IUL reserves the right to change the specifications of products at any time.

In an effort to produce useful and appropriate documentation, we appreciate your comments on these Instructions for Use. Please contact our Technical Service with any feedback.

1.6.4. Requirements for iPeak® reader 4,3" Users

The following table covers the general level of competence for the use and servicing of the iPeak® reader 4,3" .

Task	Personnel	Training and experience
Device operation	Laboratory personnel or equivalent trained operators	Trained in the procedures and handling techniques associated with the applicable assays and device operations.
Servicing	Service Specialists only	Trained and authorized by IUL to perform servicing, maintenance and repair activities.

2. SAFETY INFORMATION

Before using or servicing the iPeak® reader 4,3" , it is essential that you read this user manual carefully. Following the instructions and safety information in this user manual will ensure safe operation and maintain the system in a safe condition.

The following types of safety information appear throughout iPeak® reader 4,3" Instructions for Use.

WARNING 	The term WARNING is used to inform you about situations that could result in personal injury to you or other people. Details about these circumstances are given in a box like this one.
---	---

CAUTION 	The term CAUTION is used to inform you about situations that could result in damage to the instrument or other equipment. Details about these circumstances are given in a box like this one.
---	--

The advice given in this Instructions for Use is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

2.1. Proper Use

WARNING/ CAUTION 	Risk of personal injury and material damage Improper use of the iPeak® reader 4,3" instrument may cause personal injury or damage to the instrument. The instrument must only be operated by qualified personnel.
--	--

CAUTION 	Damage to the instrument Avoid spilling water or chemicals onto the iPeak® reader 4,3" instrument. Damage caused by water or chemical spillage will void your warranty.
---	--

In case of emergency, switch off the iPeak® reader 4,3" at the power switch and unplug the power cord from the power outlet.

2.2. Electrical Safety

If the operation of the iPeak® reader 4,3" is interrupted in any way (e.g., due to interruption of the power supply or a mechanical error), first switch off the instrument using the power switch, then disconnect the electrical cord from the power outlet. Contact Technical Service after such an incident.

WARNING 	Rechargeable Batteries The iPeak® reader 4,3" contains a rechargeable lithium-ion battery pack inside that accomplishes the IEC 62133 standard: Secondary cells and batteries containing Alkaline or other Non-Acid Electrolytes – safety requirements for portable sealed secondary cells and for batteries made from them, for use in portable applications. This WARNING is applicable to routine users as well as to technical service users. • The battery is automatically charged when the device is connected to the external power supply. • The battery is not user-replaceable. Replacement must only be performed by manufacturer-authorized service personnel. • Do not attempt to dismantle, open, replace or shred the battery. • At the end of the device lifetime, the device and its battery must be disposed of according to applicable regulations for electronic waste and batteries. • Keep batteries out of the reach of children. • Seek medical advice immediately if a battery has been swallowed. • Do not expose batteries to heat or fire. Avoid storage in direct sunlight. • Do not remove a battery from its original packaging. • Do not subject batteries to mechanical shock • In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice. • Only use Batteries provided by IUL. • After a storage period of 6 months, it is necessary to charge the batteries to obtain maximum performance.
---	--

WARNING 	Electrical hazard Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument unsafe. This must be checked after service or maintenance. Intentional interruption is prohibited. Lethal voltages inside the instrument This WARNING is applicable to routine users as well as to technical service users. • Risk of electrical shock and energy hazard. All failures should be examined by a qualified technician. Please do not remove the case of the AC adaptor by yourself! • Adaptors should be placed on a reliable surface. A drop or fall could cause damage. • Please do not place the AC adaptor in places with high moisture or near the water. • Please do not place the AC adaptor in places with high ambient temperature or near a fire source. About the maximum ambient temperature, please refer to "Appendix B: Technical Data". • Disconnect the AC adaptor from the AC power before cleaning. Do not use any liquid or aerosol cleaner. Only use a damp cloth to wipe it. In case of replacement or a loosening of the AC adaptor or the mains power cord, these must be replaced only with the AC adaptors or power cords listed in the "Ordering codes" Appendix. In case of replacement, this must be ordered to IUL.
---	---

To ensure satisfactory and safe operation of the iPeak® reader 4,3" :

- The power cord of the external power supply must be connected to a line power outlet that has a protective conductor (earth/ground).
- No other external power supply neither power cords other than the ones specified in the 'Ordering Codes' Appendix must be used. In case of replacement, this must be ordered from IUL.
- The instrument must not be operated with the cover removed.
- If you suspect any instrument damage, contact Technical Service.

If the iPeak® reader 4,3" becomes electrically unsafe, prevent other personnel from operating it, and contact Technical Service.

The instrument may be electrically unsafe if:

- The instrument or the external power supply is shown to be damaged.
- The instrument has been stored under unfavourable conditions for a prolonged period.
- A different external power supply is used other than the one provided by IUL's partner.

WARNING 	Risk of electric shock
---	--------------------------------------

	In case of replacement or loss of the external power supply or the power cord, these must be replaced only with the external power supply or power cords listed in the 'Ordering Codes' Appendix and provided by IUL's partner.
--	---

2.3. Environment

Operating conditions

WARNING 	Explosive atmosphere The iPeak® reader 4,3" is not designed for use in an explosive atmosphere.
CAUTION 	Direct sunlight Do not expose the iPeak® reader 4,3" to direct sunlight or other powerful lights during operation.
CAUTION 	High humidity or liquids Protect the iPeak® reader 4,3" from high humidity and contact with liquids.
CAUTION 	Strong electromagnetic radiation Do not expose the iPeak® reader 4,3" to strong electromagnetic radiation.
CAUTION 	Strong ultrasonic radiation Do not expose the iPeak® reader 4,3" to strong ultrasonic radiation.

2.4. Biological Safety

Use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS:

https://www.cdc.gov/labs/pdf/SF_19_308133-A_BMBL6_00-BOOK-WEB-final-3.pdf

WARNING 	Samples containing infectious agents Some samples used with the iPeak® reader 4,3" may contain infectious agents. Handle such samples in accordance with the required safety regulations. The responsible person(s) (e.g., laboratory manager) must take the necessary precautions to ensure that the workplace is safe and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents, as defined in the applicable Safety Data Sheets (SDSs) or OSHA ¹ ACGIH ² or COSHH ³ documents. Venting of fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.
WARNING 	Samples containing infectious agents Use safe laboratory procedures as outlined in publications such as <i>Biosafety in Microbiological and Biomedical Laboratories</i> , HHS:

¹ OSHA: Occupational Safety and Health Administration (United States of America).

² ACGIH: American Conference of Government Industrial Hygienists (United States of America).

³ COSHH: Control of Substances Hazardous to Health (United Kingdom).

	https://www.cdc.gov/labs/pdf/SF_19_308133-A_BMBL6_00-BOOK-WEB-final-3.pdf See: Section III—Principles of Biosafety Laboratory Practices and Techniques The most important element of containment is strict adherence to standard microbiological practices and techniques. Persons working with infectious agents or potentially infected materials must be aware of potential hazards and must be trained and proficient in the practices and techniques required for handling such material safely. The director or person in charge of the laboratory is responsible for providing or arranging the appropriate training of personnel.
WARNING 	Samples containing infectious agents To avoid hazards in case of cassette break use Primary Barriers and Personal Protective Equipment: Use safe laboratory procedures as outlined in publications such as <i>Biosafety in Microbiological and Biomedical Laboratories</i>, HHS: https://www.cdc.gov/labs/pdf/SF_19_308133-A_BMBL6_00-BOOK-WEB-final-3.pdf See: Section III—Principles of Biosafety Safety Equipment (Primary Barriers and Personal Protective Equipment) Safety equipment includes BSCs, enclosed containers, and other engineering controls designed to remove or minimize exposures to hazardous biological materials.

2.5. Chemicals

WARNING 	Hazardous chemicals Some chemicals used with the iPeak® reader 4,3" may be hazardous. Always wear safety glasses, gloves, and a lab coat. The responsible person(s) (e.g., laboratory manager) must take the necessary precautions to ensure that the workplace is safe and that the instrument operators are suitably trained and not exposed to hazardous levels of toxic substances (chemical or biological), as defined in the applicable Safety Data Sheets (SDSs) or OSHA¹ ACGIH² or COSHH³ documents. Venting of fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.
---	--

¹ OSHA: Occupational Safety and Health Administration (United States of America).

² ACGIH: American Conference of Government Industrial Hygienists (United States of America).

³ COSHH: Control of Substances Hazardous to Health (United Kingdom).

2.6. Maintenance Safety

CAUTION 	Damage to the instrument Do not use spray bottles containing acetone or disinfectant to clean the surface of the iPeak® reader 4,3" instrument. Do not use products containing acetone or other corrosive solvents to clean the iPeak® reader 4,3" instrument.
WARNING 	Risk of electric shock Do not open the panels on the instruments. Risk of personal injury and material damage Only perform maintenance that is specifically described in these Instructions for Use.

2.7. Waste Disposal

Used consumables, such as cassettes, may contain hazardous chemicals or infectious agents. Collect and dispose of them following local safety regulations.

All packaging waste from cassettes must be collected and disposed of properly following local environmental regulations.

For disposal of waste electrical and electronic equipment (WEEE) see the corresponding appendix.

2.8. Glossary of symbols

The symbols shown in the table below **may or may not appear** on your device. The presence of specific symbols depends on applicable regulatory requirements and the country or region where the product is marketed.

Symbol	Location	Description
	Type plate on the bottom of the instrument and the external box label	CE Mark, Declaration of Conformity, see the corresponding appendix.
	Type plate on the bottom of the instrument and the external box label	UKCA Mark, Declaration of Conformity, see the corresponding appendix.
	Type plate on the bottom of the instrument and the external box label	QPS certification mark: The product has been tested and found in compliance with the relevant product electrical safety standards of both the U.S. and Canada.
	Type plate on the bottom of the instrument and the external box label	Medical device manufacturer.

Symbol	Location	Description
	Type plate on the bottom of the instrument and the external box label	The device is intended to be used as an in vitro medical device.
	Type plate on the bottom of the instrument and the external box label	Date of manufacture.
	Type plate on the bottom of the instrument and the external box label	Catalog Number.
	Type plate on the bottom of the instrument and the external box label	Serial Number.
	Label on the side of the external box	Indicates a carrier that contains unique device identifier information.
	Label on the side of the external box	Temperature limits, see 'Technical Data' Appendix.
	Label on the side of the external box	Humidity limits, see 'Technical Data' Appendix.
	Only in the cassettes that have a biological hazard (information on the cassettes' instructions for use)	Biological hazard symbol.
	Type plate on the bottom of the instrument and the external box label	Consult instructions for use.
	Type plate on the bottom of the instrument	Caution is necessary when operating the device.
	Type plate on the bottom of the instrument and the external box label	FCC compliant, Declaration of Conformity, see the corresponding appendix.
	Label on the side of the external box	All components of the instrument are recyclable.
	On the bottom of the drawer	Identify the material the drawer is made of (Polyamide).

Symbol	Location	Description
	Type plate on the bottom of the instrument and the external box label	Waste Electrical and Electronic Equipment (WEEE), see the corresponding appendix.

2.9. Cybersecurity

The LFR Software is an essential part of the LFR and it is not distributed separately.

All the cybersecurity controls are integrated into the SW storage, including SW integrity verification, firewall, and others.

The LFR is designed to operate in the intended environment of use defined by IUL's partner. The LFR's partner is the legal entity that defines the setup for a specific intended environment of use, including the LFR's setup and operation process.

The device is intended to operate within a controlled and secured network environment managed according to the user organization's IT and cybersecurity procedures.

For this purpose, the LFR provides a privileged user (role Administrator) which is the only user that could configure 'Cassette' and 'Preferences' tabs.

There is another kind of privileged user (Manager role) that can manage users, printers, and network setups. Only the Administrator and Manager have access to this feature.

The LFR starts with a Login Screen, to introduce the user and password before operating.

Strong passwords and restricted administrator access are recommended.

There is an audit trail integrated into the SW to store performed actions involving medical data.


Software updates provided by the manufacturer or its authorized partner should be installed according to the instructions provided.

NOTE: Contact Technical Support to request technical specifications of Remote Export and Audit Trail functionalities.

3. GENERAL DESCRIPTION


3.1. Device Overview

Front View:

- 1. Touchscreen and display
- 2. Barcode reader
- 3. Drawer insertion slot
- 4. Power ON indicator (white LEDs)
- 5. Power button
- 6. Battery in charge red LED

Rear View:

- 7. External power supply socket
- 8. Ethernet socket
- 9. USB host port (only accessible for slim USB memory sticks)
- 10. Software Wipe button

CAUTION

For proper insertion of a slim USB memory stick, its body section near the connector must be 8,5x16mm or less.

3.2. Accessories of the device

3.2.1. Drawer

The drawer model must fit with the corresponding Cassette model for proper use of the LFR. See the 'Ordering codes' appendix to find all available models.

3.2.2. External Power Supply

Connect the external power supply to mains, and the LFR to the external power supply to charge it. The red LED above the power supply socket should be lit when the power supply is properly connected.

WARNING

Risk of electric shock

In case of replacement or loss of the external power supply, this must be replaced only with one of the power cords listed on the 'Ordering codes' appendix delivered by IUL.

3.2.3. Power cord

The LFR is equipped with a power cord with a plug suitable for the destination country. See the 'Ordering codes' appendix for details.

WARNING 	Risk of electric shock In case of replacement or loss of the power cord, this must be replaced <u>only</u> with one of the power cords listed on the 'Ordering codes' appendix delivered by IUL.
---	--

3.2.4. Adjustment drawer

The adjustment drawer is an accessory to the iPeak® reader 4,3" and may be used to readjust the camera parameters. It is intended solely for corrective maintenance in case of camera parameter deviation. It is not required for routine operation, preventive maintenance, or daily operation.

To perform the adjustment procedure, follow the instructions displayed on the screen step by step. The automatic adjustment process must not be interrupted once initiated.

The adjustment drawer **bears its own expiration date**. Use beyond its expiration date is not permitted and will prevent execution of the adjustment procedure.

Expiration of the adjustment drawer does not limit the service life of the iPeak® reader 4,3" .

The adjustment drawer shall be stored and handled in accordance with the storage conditions indicated on its label. Failure to respect the specified storage conditions may affect its performance and compromise the adjustment procedure.

WARNING 	During the adjustment process, the battery level should be above 70%. If the battery is too low, connect the LFR to a charger. Keep the adjustment drawer surface clean and free from dust. Use only a rubber air blower or equivalent non-contact cleaning method.
---	--

CAUTION 	Internal temperature To use the adjustment functionality, the internal temperature of the iPeak® reader 4,3" should be stable. To stabilize the internal temperature, the device should be fully charged and switched on for at least 30 minutes. The ambient temperature for the adjustment process should be $22.5 \pm 2^{\circ}\text{C}$.
---	---

4. INSTALLATION

4.1. Unpack the device

The packaging of the iPeak® reader 4,3" can be stored for reuse.

CAUTION

Before switching on the reader for the first time or after a long period of inactivity, it is essential to plug it into the mains supply. Failure to do so may cause permanent damage to the unit.

4.2. Site requirements

Place the iPeak® reader 4,3" device on a stable surface, far from powerful lights and near an earthed/grounded electrical outlet.

It is recommended to connect the iPeak® reader 4,3" device to a power supply when used continuously.

CAUTION

The device shall be used under laboratory standard conditions, in an ambient temperature between 20°C and 25°C.

4.3. Power cable connection

The socket used to connect the power cable is located at the back of the iPeak® reader 4,3" device.

When the iPeak® reader 4,3" is not in use for a long period of time, it is recommended to disconnect the power cable.

4.4. Cassette Requirements

Before using the iPeak® reader 4,3" device make sure that the cassette has not exceeded its expiration date.

WARNING

Do not use expired cassettes. Review the expiration date of the cassette and do not use it if it is expired.

5. OPERATING PROCEDURES

This section outlines the functionalities of the LFR. Note that the appearance of screenshots may vary depending on your device model.

CAUTION


Before switching on the reader for the first time or after a 6- month period of inactivity, it is essential to plug it into the mains supply. Failure to do so may cause permanent damage to the unit.

5.1. Login screen

To start the LFR, press the power button and wait until the Login Screen is shown:

To operate the LFR, properly fill in the 'Username' and 'Password' fields and then press the 'Ok' button. To edit each field, click on the field to access the virtual keyboard. Type the correct username and/or password and press the 'Enter' button. Special characters cannot be used.

NOTE: The default operator username is "operator" and its default password is "pwd".

CAUTION

Default password MUST be changed to preserve the correct use of the LFR.


See the "User Management" section for explanations on how to change passwords. Retain the new password in a secure location.

NOTE: If you forget or lose your password, contact your LFR's Manager.

To end the user's session, select 'Logout' from the Tools Menu or switch off the LFR using the power button.

5.2. Home screen overview



The Home screen displays the current Cassette Settings, date, time, battery status, 'New scan' button, and 'System menu' button.

Main screen:

	Starts a new Cassette Processing.
	Enters the System Menu.


Top bar:

	Shows the Side Menu.
	Returns to Home screen.
	Shows the information of the current Cassette Settings available to be used.
	Shows the drop-down Tools Menu.
	Shows the date, time, and battery status.
	Airplane icon. Only visible when airplane mode is enabled.

5.3. Cassette Processing

To process a cassette, follow the next steps:

1. **New Scan:** To start a new cassette scan, press the 'NEW SCAN' button on the Home screen.

2. **Cassette Settings:** Skip this step if the Cassette Settings displayed on the Top Bar is the Cassette Settings to be used in the next scan. This step may be required in certain configurations.

To set the proper Cassette Settings that will be used in the next scan, one of the three following options can be performed:

- Option A: Load a new Cassette Settings through the barcode reader.
- Option B: Manually select an available Cassette Settings from the list.
- Option C: Automatically select an available Cassette Settings from the list.


Below are the detailed steps for each of the three options:

Option A: load a new Cassette Settings through the barcode reader.

To load a new Cassette Settings, press the 'Cassette config' button. Depending on the LFR configuration, the screen is displayed in two different ways:

Take the Data Matrix of the correct Cassette Settings, which is supplied along with the cassette box:

Place it in front of the LFR's barcode reader approximately 5 cm away. Press the 'Activate barcode read' button. The barcode reader starts its reading process using a light. This light disappears when a barcode has been read or after 10 seconds of trying.


NOTE: If a barcode cannot be read, check and clean the window of the barcode reader.

If the barcode is read but does not fit with the expected features, an error message is shown on the screen. Check the scanned barcode and press any part of the display outside the dialog box to remove the error message and return to the previous screen.

If the barcode is read correctly, the new Cassette Settings data is shown on the Top Bar.

Option B: Manually select an available Cassette Settings from the list.

If the Cassette Settings to be used was previously uploaded, it is not required to upload the Cassette Settings again. Go to the Side Menu, press 'Cassette settings', and the list of Cassette Settings is displayed.

To load the desired Cassette Settings, press the corresponding 'Select as current' button.

Option C: Automatically select an available Cassette Settings from the list.

If the system is ready to work with "Barcode Autodetection" and the cassette has the proper embedded barcode, it is possible to directly launch the scan and the LFR will use the Cassette Settings of the embedded barcode as the current Cassette Settings.

Pre-requisites:

- The inserted cassette has the proper embedded barcode.
- The corresponding Cassette Settings is available on the list.
- There is a valid (not expired) Cassette Settings displayed on the Top Bar.

3. **Confirmation:** Once the proper Cassette Settings is set, press the 'Confirm cassette' button to check that the current Cassette Settings match the cassette to be processed. Depending on the configuration, this step may be disregarded, and the button may not be present.

Place the small barcode (Small-Data Matrix) in front of the barcode reader at 5 cm of distance approximately.


Press the 'Activate barcode read' button. The barcode reader starts reading using the light. This light disappears when a barcode has been read or after it has tried reading for 10 seconds.

If the Small-Data Matrix does not correspond to the Cassette Settings currently loaded in the LFR, an error message is displayed. In this case, use a cassette of the proper type or load another Cassette Settings corresponding to the current cassette.

WARNING 	It is the operator's responsibility to use the Small-Data Matrix (cassette's confirmation barcode) from the current cassette.
---	---

4. **Sample ID:** After reading a Data Matrix, the LFR returns to the scanning screen and the Sample ID may be entered. Depending on the LFR configuration, this step may be mandatory. There are two options to enter the Sample ID:
- Option A: Manually type the Sample ID
 - Option B: Read the sample ID from a barcode

Depending on the LFR configuration, the screen appears in 3 different ways:


Below are the detailed steps for each of the two options:

Option A: Manually type the Sample ID.


Click on the 'Enter sample ID' field and a virtual keyboard will appear. After typing the sample ID with the virtual keyboard, press the 'Enter button to proceed.

NOTE: Special characters cannot be typed in and the maximum length of the text is 30 characters.



Option B: Read the sample ID from a barcode.

With this option, special characters are allowed and the maximum length is 1024 characters (if the number of characters exceed 1024, the subsequent characters will be omitted). To load the sample ID from a barcode, press the 'Activate barcode read' button and place the barcode in front of the light.

In both options and depending on the LFR configuration, if the Sample ID entered is already present in the stored results of the LFR, a message is displayed, and the scan button may be disabled.



Go the "Scan Results" section to see how Sample IDs are displayed in different scenarios.

- 5. Drawer:** Insert the cassette into the corresponding drawer. Ensure that the code in the drawer bottom (see green circle) matches with the drawer code displayed along with the cassette settings at the Top Bar. Ensure that the cassette and the drawer surfaces are levelled. Push the drawer inside the LFR.

WARNING	Before inserting a cassette:
	Ensure that the four (4) white dots (fiducials) on top of the drawer are present, undamaged, clean (see the "Maintenance" section), and round (i.e., in their original factory conditions).
	Ensure that strips embedded in the cassette are free of impurities such as fibres or dust.
	Failing to ensure any of the above may lead to an incorrect result.

A warning is shown when the drawer is not properly closed. Depending on the LFR configuration, this screen appears in two different ways:

The drawer must be inserted correctly and completely for the error message to disappear.

- 6. Launch:** The scan buttons ('Scan' and 'Timed scan') remain in grey until some conditions are reached, depending on the LFR configuration:
 - A valid Cassette Settings is loaded
 - The Cassette Settings is confirmed (if required)

- A Sample ID is set (if required)
- A drawer is placed inside the LFR.

Once the required conditions are met, the scan buttons go from grey to colored. Two options are available to process a cassette: an immediate scan ('Scan' button) or a timed scan ('Timed Scan' button).

- Option A: press the 'Scan' button to launch a scan immediately
- Option B: press 'Timed Scan' button to delay the scan launch as prescribed.

Depending on the LFR configuration, this screen will appear in two different ways:


Below are the detailed steps for each of the two options:

Option A: Press the 'Scan' button to launch a scan immediately.

If the prescribed time after sample inoculation has passed, press the 'Scan' button to scan the cassette immediately.

Option B: Press the 'Timed Scan' button to delay the scan launch as prescribed.

If the sample has just been inoculated, press the 'Timed Scan' button to defer the scan to the prescribed waiting time. After pressing the button, a countdown is shown until the scan starts. Press the 'Cancel' button to stop the process.

- 7. Scan result:** The results of the scan are displayed when the process ends. Scroll up and down on the screen and press the different tabs to see all the data. See 'Scan Results' section for more information about Results Screen and the functionalities of each button.

STRIP_01	IMAGES	COMMENTS	INFO
test 1: 52.8 ng/ml test 2: Positive (1234) Control Line: Present (1403)			
PRINT TICKET EXPORT PDF PRINT PDF NEW SCAN			

To process another cassette, press the 'New scan' button on the bottom-right corner.

WARNING 	The reader will provide reliable results if the tests are performed according to the IUL's partner specifications. The control and test lines are differentiated from the background and the cassette does not present artifacts or impurities such as fibre or dust.
--------------------	---

5.4. System Menu

Set the LFR configuration through this menu. Different tabs are available, depending on the user's role.

5.4.1. System

NOTE: Depending on the LFR configuration some items are not displayed.

Language	English ▼
Format USB	FORMAT
Export Logs to USB	EXPORT
Set Time Zone	SET TIME ZONE
Remote Support	CONNECT
Adjust Capture Parameters	ADJUST
Activate Line Quality	READ QR

Language:	Use this to select a language among those installed in the LFR.
Inactivity Shutdown Timeout (min):	Time in minutes from the last interaction with the LFR until the system powers down itself to save battery. This parameter may or may not be present, depending on the LFR's configuration. The default value is 30 minutes. Contact your Administrator to change its maximum value or disable the functionality.
Format USB:	Press the "Format" button to format in FAT32 any USB detected. This results in the removal of previous content. The USB can then be used to export CSV and PDF results, write logs, and update the LFR software.
Export Logs to USB:	Press the "Export" button to write internal data to a slim USB memory stick plugged into the USB host port. Important for the Technical Service as an aid for diagnosis.
Set Date/Time:	Press the button to define the current date and time. This parameter has configurable access. NOTE: The LFR cannot work with dates before the manufacturing date.
Set Time Zone:	Press the button to define the current time zone.

CAUTION 	After changing the Time Zone, re-start the LFR.
---	---

Remote Support: Press the 'Connect' button whenever technical support is needed.

CAUTION 	To use the 'Connect' button. Internet connection is needed, specifically bidirectional connectivity through port 22 with the Technical Support.
---	---

Adjust capture parameters:

Press this button to readjust the camera parameters. To perform the adjustment, first charge the LFR's battery to 100%, and leave it connected to the main power source. To use the adjustment functionality, the internal temperature of the device should be stable. To stabilize the internal temperature, the device should be fully charged and switched on for at least 30 minutes. The ambient temperature for the adjustment process should be $22.5 \pm 2^{\circ}\text{C}$. Insert the adjustment drawer (Cat. 900013732) and press 'Adjust' to start the process.

WARNING 	The adjustment drawer must be in good condition, otherwise, the process could misadjust the LFR. Avoid direct illumination on the drawer slot.
---	--

Activate Line Quality:

Press the button to activate the Line Quality functionality by scanning a code with the barcode reader.

NOTE: This code is serial number-specific. This code is provided by IUL.

IMPORTANT: The Line Quality functionality is an additional option to the standard software. Contact sales@iul-inst.com for more information.

5.4.2. Cassette (Administrator only)

Only the Administrator user can access this screen. Parameters related to Cassette Settings are displayed and can be edited:

Cassette Settings Version
0

0-1000

SET

Keep Cassette Settings when Power Down

☒

Maximum Cassette Settings

10

Require Confirmation

☒ None
☐ Small Barcode
☐ Barcode Autodetection
☐ Configuration Barcode

Cassette Settings Version:

This parameter is used to accept or block some Cassette Settings depending on its configuration. This is an internal parameter of the LFR, but the Cassette Settings also has a parameter called Cassette Settings Version. When the Cassette Settings Version of the LFR is set to "0", the LFR accepts **all** Cassette Settings. When the Cassette Settings Version of the LFR is set to a value other than "0" for example "3", the LFR accepts Cassette Settings Versions **only** with the numbers "3" or "0" and **blocks all others**. The Cassette Settings Version of the Cassette Settings is set when exporting the Cassette Settings from the Development Software Application. Below are some examples of different combinations:

Value of the Cassette Settings Version in iPeak® reader 4,3"	0	0	7	7	3
Value of the Cassette Settings Version in Cassette Settings:	0	3	0	3	3
Does the LFR accept the Cassette Settings?	YES	YES	YES	NO	YES

Keep Cassette Settings when Power Down:

If ON, the current Cassette Settings is kept when the LFR is switched off; otherwise, the current Cassette Settings is lost when the LFR is switched off.

Maximum Cassette Settings:

Number of Cassette Settings that can be stored in the LFR and available to use (without re-reading its corresponding Data Matrix). The maximum value to set is 20 and the minimum is 1.

Require Confirmation:

Four options are available: 'None', 'Small Barcode', 'Barcode Autodetection', and 'Configuration Barcode'. Below are details of each option:

- *None*: It is possible to process samples with current Cassette Settings without restrictions.
- *Small Barcode*: The corresponding confirmation barcode (Small-Data Matrix) must be read before each scan, to confirm that the current cassette matches the current Cassette Settings. Use the frontal barcode reader to do so (see 'Device Overview' section).
- *Barcode Autodetection*: After launching a scan (with 'Scan' or 'Timed scan' buttons), the LFR will try to find the Cassette Settings that match one on the Cassette Settings list by reading the barcode embedded on the cassette. Using this option, you can also read the Small-Data Matrix with the barcode reader to (see 'Device Overview' section).

CAUTION 	Default password MUST be changed to preserve the correct use of the LFR.
---	---

- *Configuration Barcode*: It is necessary to load the proper Cassette Settings before each scan.

5.4.3. Preferences (Administrator only)

Only the Administrator user can access this screen. Some parameters are displayed and could be edited:

Sample ID Checking	☐ No Check ☐ Required ☐ Warn if Duplicated ☒ Block if Duplicated
Maximum Results	200 SET
Check Database is Full	☒
Current Scan Counter: 2	RESET
Inactivity Shutdown Enabled	☒
Inactivity Shutdown Time Max. (min)	30
Set Date/Time Enabled	☒
Remote Export Enabled	☒
HL7 CR: 0D	☒
HL7 Export Message Type	MDM_T02
SFTP Export Enabled	☒
TCP/IP Export Enabled	☒
TCP/IP ACK Enabled	☒
Reset Cassette Settings Enabled	☒
Show Peak Area to Operator	☒
Show Profile Image to Operator	☒
Show Expiration Date to Operator	☒

Sample ID Checking:

Four ID checking levels are available: 'No Check', 'Required', 'Warn if Duplicate', and 'Block if Duplicate'.

- *No Check*: There are no restrictions or warnings.
- *Required*: A value is mandatory for the scan to proceed.
- *Warn if Duplicated*: A value is mandatory for the scan to proceed, and a warning is displayed if the current value is already in the LFR's stored results.
- *Block if Duplicated*: A value is mandatory for the scan to proceed, and the 'Scan' and 'Timed scan' buttons are disabled if the current value is already in the LFR's stored results

**Maximum
Results
Database:**

Set the maximum number of results that are stored in the LFR (up to 200). When the limit is exceeded, the oldest data is overwritten.

WARNING 	To avoid losing older results, before decreasing the Max Results Database, an Export Results procedure must be performed (see 'Scan Results' section).
---	--

**Check Database
is Full:**

When the database reaches 80% of its capacity, a warning message is shown when the LFR is powered on and when a new Cassette Processing is started. When the message shows 100%, the oldest result is overwritten when a new scan is performed. This option can be disabled if storing results is not mandatory/important.

**Current Scan
Counter:**

Displays the number of scans that have been performed since the last reset. Press the 'Reset' button to reset this counter.

**Inactivity
Shutdown
Enabled:**

If ON, the LFR switches off automatically when the prescribed time is surpassed without any activity detected.

**Inactivity
Shutdown Time
Max. (min):**

Maximum time in minutes that the operator can set if the inactivity shutdown is enabled.

**Set Data/Time
Enabled:**

If ON, it is possible to set Date/Time in the 'System' tab.

**Remote Export
Enabled:**

Enables the remote data export functionality. If enabled, an external server can send a series of requests to the LFR that will be answered.

HL7 CR: 0D:

Defines the segment separator used in HL7 message export:

- *Enabled:* The system uses the non-printable character for carriage return as the HL7 segment separator. For 'Remote Export' functionality, the message data is base64 encoded to ensure compatibility with web communication protocols. For 'SFTP Export' and 'TCP/IP Export' functionalities, a file is sent whose data is not encoded.
- *Disabled:* The system uses the four-characters string "<cr>" as the HL7 segment separator.

**HL7 Export
Message Type:**

Two options are available:

- *ORU_R01:* ORU_R01 result message structure
- *MDM_T02:* MDM_T02 result message structure

**SFTP Export
Enabled:**

If ON, after each scan the SFTP export service tries to send the result to its configured SFTP server. If an error is found a pop-up dialog shows in its description.

When SFTP export is enabled, the Scan Results list displays an SFTP status for each scan. In addition, if there is a configured SFTP server, a 'Retry SFTP' button is available in case of an error.

TCP/IP Export Enabled:	If ON, after each scan the TCP/IP export service tries to send the result to its configured TCP/IP server. If an error is found a pop-up dialog shows its description. When SFTP export is enabled, the Scan Results list displays an TCP/IP status for each scan. In addition, if there is a configured TCP/IP server, a 'Retry TCP/IP' button is available in case of an error. NOTE: Contact Technical Service to request technical specifications of results export functionalities.
TCP/IP ACK Enabled:	If ON the acknowledgement reply shall be used in the TCP/IP communication.
Reset Cassette Settings Enabled:	If ON, users with the Operator role can restore the Cassette Settings to their factory default. If no factory cassette settings are provided in the reader, the list will be cleared/emptied.
Show Peak Area to Operator:	If ON, users with the Operator role will see the peak area displayed in brackets for qualitative results. When set to OFF and the logged user is an Operator, the peak area will be hidden when the result is displayed or exported (ticket, PDF report, and CSV file).
Show Profile Image to Operator:	If ON, users with the Operator role will be able to view the profile image in the Image tab of the results. When set to OFF and the logged user is an Operator, the profile image will be hidden when the result is displayed or exported PDF reports.
Show Expiration Date to Operator:	If ON, users with the Operator role will see the expiration date of the cassette settings introduced in the reader. When set to OFF and the logged user is an Operator, the expiration date will be hidden when the result is displayed or exported (ticket, PDF report, and CSV file).

5.4.4. Printer (Administrator or Manager only)

Only users with the roles of Administrator or Manager can access this screen.

Desktop printer	HP_DeskJet_2640 ▾
Baud Rate	9600 ▾
Parity	None ▾
Data Bits	8 ▾
Stop Bits	1 ▾

A desktop printer model could be selected, and communication parameters for a ticket printer are displayed and could be edited. Both the desktop printer and ticket printer must be connected to the device's USB host port (see 'Device overview' section) for instructions on how to use them.

Desktop printer: Allows the selection of the printer for printing reports in PDF format.

To use a ticket printer, the following values corresponding to the configuration of the ticket printer to must be entered:

- Baud rate:** Method for detecting errors in the transmission of serial communication with the ticket printer.
- Parity:** Method for detecting errors in the transmission of serial communication with the ticket printer.
- Data bits:** Number of bits per character of the serial communication with the ticket printer.
- Stop bits;** Number of stop bits sent at the end of each character of the serial communication with the ticket printer.

The ticket printer MUST support ESC/POS control language, and there are two options:

- Option A: Ticket printer with serial port, and an external USB to the serial port adaptor.
- Option B: Ticket printer with USB port, internally emulating serial port.

5.4.5. Network (Administrator or Manager only)

Only users with the roles of Administrator or Manager can edit the parameters of this screen.

SYSTEM	CASSETTE	PREFERENCES	PRINTER	NETWORK	VERSIONS
Airplane Mode		☐			
Ethernet		SET			
Status		On			
MAC		XX:XX:XX:XX:XX:XX			
Mode		DHCP Client			
IP		XXX.XXX.X.XXX			
Mask		XXX.XXX.X.XXX			
Broadcast		XXX.XXX.X.XXX			
Gateway		XXX.XXX.X.XXX			
DNS		XXX.XXX.X.XXX			
WLAN		SET			
Status		On			
MAC		XX:XX:XX:XX:XX:XX			
Mode		Access Point			
SSID		lfr_xxxx			
IP		XXX.XXX.X.XXX			
SFTP		SET			
Not defined					
TCP/IP		SET			
Not defined					
Forward WLAN to Ethernet		☐			

NOTE: The labels shown in this figure are for illustrative purposes only. Values are displayed in a generic format (e.g., XXX.XXX).

Users with the Operator role can only enable or disable the Airplane Mode. Network communications parameters are displayed:

Airplane Mode: If ON, Wi-Fi connections are disabled, and an airplane icon is shown on the Top Bar.

Ethernet: Displays the current Ethernet configuration. Press the ‘Set’ button to configure it. There are two modes:

- *DHCP Client mode:* The network server assigns network parameters to the LFR. This is the default mode.
- *Static mode:* The proper connection parameters must be set on the LFR (specifically: IP, Mask, and Gateway).

Mode

Static

IP

Mask

Gateway

CLOSE

SAVE

Mode

DHCP Client

CLOSE

SAVE

WLAN: Displays the current Wi-Fi network’s configuration. Press the ‘Set’ button to configure it. There are two modes:

- *Access Point mode:* This mode allows the LFR to function as an access point. This is the default mode. It requires defining an SSID and a password, as well as selecting a subnet from those available.
- *DHCP Client mode:* This mode allows the LFR to connect to an external Wi-Fi network by entering the SSID and password.

Mode

Access Point

SSID

lfr_2037

Password

Subnet

192.168.138.0

CLOSE

SAVE

Mode

DHCP Client

SSID

lfr_2037

Password

CLOSE

SAVE

NOTE: The WLAN SSID field must contain between 1 and 32 characters and the WLAN Password field must contain between 8 and 63 characters (IEEE 802.11i-2004 compliant). To edit the WLAN configuration, the Airplane Mode must be disabled.

SFTP: Displays the current SFTP server configuration. Press the ‘Set’ button to configure it.

Server

Port

22

HL7 path


PDF path

Username

Password

SAVE

TCP/IP: Displays the current TCP/IP server configuration. Press the ‘Set’ button to configure it.

NOTE: When setting up some of these servers, if the LFR has the corresponding functionality activated, it will attempt to export all results that have not been previously exported.

Forward WLAN to Ethernet: If enabled, other LFRs connected to its Access Point will send data to server through its Ethernet connection.

DISCLAIMER:
The use of the LFR as a forwarding hub for other devices may be subject to restrictions derived from the licensing conditions of the end user’s LIS system. IUL does not guarantee the legal or technical compatibility of this functionality with third-party laboratory information systems. Both the distributor/integrator of the solution (OEM partner) and the end user assume full responsibility for validating this integration with the LIS provider and for obtaining any additional licenses that may be required. Under no circumstances shall IUL be held liable for any improper or unauthorized use of the system in this context.

5.4.6. Versions

Displays a list of information related to the LFR and its software.

Model	Device name
Serial Number	1000xx000/xxxxx
Current SW Version	x.x
Hardware Version	xxxxxx
Product Image Version	abcde-fghij-x.x-s
Wipe Clean SW Version	x.x
Factory SW Version	x.x


NOTE: The labels shown in this figure are for illustrative purposes only. Values are displayed in a generic format (e.g., XXX.XXX).

Model: Commercial name of the LFR.

Serial Number:	Unique number for each LFR. This can be found at the label at the bottom of the LFR.
Current SW Version:	Software version of the LFR.
Hardware Version:	Hardware version of the LFR.
Product Image Version:	Name of image file initially saved in current SD.
Wipe clean SW Version:	Software version if a wipe clean is performed on the LFR.
Factory SW Version:	Software version used when manufacturing.

5.5. Side Menu

By pressing the menu symbol on the left of the Top Bar, a drop-side menu is displayed as below, with the following options:

5.5.1. Main Menu

Returns to the Home screen.

5.5.2. Scan Results

Access the list of stored scan results. The list is in the order of scans, with the most recent scan on top, A maximum of 5 results are shown on each page. Scroll up and down on the touchscreen to browse the list of the results and press on each one to explore the result. At the bottom of each page, there are backward and forward buttons to change pages. On the upper-right corner, there is an 'Export Results' button to export all stored results.

8411772557210

Date: 2024-11-12 17:19 - Type: Colors Method

rf

Date: 2024-11-12 16:11 - Type: Colors Method

8411772557210

Date: 2024-11-12 13:51 - Type: Colors Method

esf

Date: 2024-11-12 12:32 - Type: Colors Method

ffyrf

Date: 2024-11-12 12:03 - Type: Colors Method

16-20 of 22

WARNING

By default, a maximum of 200 results are stored in the LFR. If you do not want to lose your old results before archiving this amount, you must export the results.

For each scan in the list, the following information is displayed:


- Sample ID reference (only first 30 characters are displayed)
- Date and time of the scan
- Cassette Settings type
- SFTP state (only when SFTP is enabled)
- TCP/IP state (only when TCP/IP is enabled)
- The state of the remote export services could have the following values:
 - *Pending*: Data sending is ongoing or there is no server configured.
 - *Success*: Scan data has been properly sent.
 - *Error*: Scan data has not been properly sent. A RETRY button is shown next to the item.

NOTE: After pressing the 'Retry' button, refresh the screen to ensure that the state of the remote export service has been properly updated.

NOTE: The results already successfully sent cannot be sent again.

1. RESULT VIEW

Select a scan result to explore its information, print it as a ticket or as a report, export it as a report or delete the whole scan result. The delete option is only for users with the Administrator or Manager roles.

Available tabs:

Strip(s):

This tab shows the names of the strips defined by the Cassette Settings. There is one 'Strip' tab for each strip defined in the Cassette Settings used for the scan. If the control line is present, result values for each test line and control line are shown. If the control line is negative, only the control line result is displayed.

NOTE: If the Line Quality functionality is enabled, the following warning may appear next to the control or test line: "Warning: (4) Invalid test: Possible artifact detected on <line_name> line". This warning may appear if one or more parameters from Line Quality checking are outside the expected range.

Contact Technical Support for more information and help with this functionality.

Images:

This tab shows processed pixels of each test line with the peak area curve profile over it. It contains one image for each strip defined on the Cassette Settings used for this scan.

WARNING 	In the Images tab, ensure that the images are aligned, focused, and not distorted. These may lead to an incorrect result.
---	--

Comments:

This tab is an open space to add comments to the scan result.

Info:

This tab displays data related to the scan process.

Available buttons:

Print ticket:

Press this button to obtain a ticket/label containing the identification data and the summary results of the scan. This can only be done if a ticket printer is connected and properly configured. See 'Printer' section for the setup and configuration of the ticket printer (Administrator or Manager role only).

```

2025-08-27 12:28:11 (CEST)
ID: sample ABC
Qualitative Method (384)
Lot: 03A9
Expiration: 2027-06-30 02:00:00 (CEST)
Reader: 100033000/99999
Username: operator
Comment:

>strip1
Control: Present (812)
Test: Positive (583)

```


Export PDF: Press this button to export the results data into a PDF report file. A slim USB memory stick (previously formatted as described in the 'System' section) must be connected to the LFR's USB host port.

NOTE: The PDF report filename contains a maximum of 228 first characters from Sample ID followed by a timestamp. Special characters are removed from the file name; spaces are replaced by dashes and uppercases are replaced by lowercases.

Print PDF: Press this button to print the scan result report. A desktop printer must be connected and properly configured. See the 'Printer' section for the setup and configuration of the desktop printer (Administrator or Manager role only).


New scan: Press this button to start a new Cassette Processing.

Delete: Press this button to completely remove the scan result from the database. This button is only available for users with the role of Administrator or Manager. When the trash icon is pressed, a pop-up dialog asks for deleting confirmation.

2. Export Results Button

After pressing the 'Export Results' button from the Scan Results screen, a pop-up dialog is displayed to select from the different export options. Before pressing the 'Ok' button a slim USB memory stick (previously formatted as described in the 'System' section) should be plugged into the USB host port.

The export process creates a folder in the root of the USB, named as current timestamp (YYYYMMDD_hhmmss) to store exported file(s). Scan Results data are transferred as a CSV file, which can be later opened with a spreadsheet application on a computer. Depending on the selected options the process performs some extra actions or not:

Include Images on Export:	If enabled, captured images from each scan are also exported in PNG format. Image file names are composed of the Sample ID (at most the first 30 characters), an underscore, and the uuid of the scan.
Purge after Export:	If enabled, the LFR's results database is removed after successful exportation.

About the exported CSV file:

- The file name starts with "Results_of_" and ends with the LFR serial number.
- The file extension is .CSV.
- The field separator/delimiter is the character TAB, and numeric decimal separator is always the character FULL STOP ".". This must be taken into consideration when importing the .CSV file into a spreadsheet application.
- The first row contains the headers/titles of the fields.
- For each scan result, there are as many rows as lines (Control+Test) in all the strips of the cassette. For instance, in the case of one cassette with 2 strips, and each strip with 1 Control line and 2 Test lines, the result is 6 rows in the .CSV.
- In the case of a strip where the Control line is not found, no rows corresponding to Test lines are present (as they cannot be evaluated if the Control line is missing).

Parameters available in the results (for each processed line):

Reference:	Sample ID (maximum length 1024 characters)
Uuid:	Unique identifier of the scan (also included in the captured image name)
TimeStamp:	UTC time of scan
Comments:	Comments added after the scan
Cassette type:	Name of Cassette Settings
Cassette code:	Code of Cassette Settings
Batch ID:	Batch ID of the cassette
Due date:	Expiring date of cassette batch
Username:	Name of the operator that performed the scan
NormError:	Processing parameter
LineError:	Processing parameter
Line warnings:	Processing parameter
Strip name:	Name of the strip from a cassette
Line name:	Name of the line from a strip

Result qualifier:	Qualification result of processed line
Quantification:	Quantification result of processed line (decimal separator ".")
Quantification units:	Units of the quantification result
Peak area:	Measured intensity of the line
Base line:	Processing parameter
Peak height:	Processing parameter
Peak position (mm):	Processing parameter
Drawer intensity:	Processing parameter
Wb_factors:	Processing parameters
Applied thresholds:	Cassette Settings parameter
Fiducials:	Processing parameter.


The following parameters are only available if the Line Quality functionality is enabled:

Linearity_horizontal:	Line quality metrics
Linearity_vertical:	Line quality metrics
Peak_symmetry:	Line quality metrics

NOTE: Depending on the LFR configuration, some fields will not be present if you have logged in with Manager or Operator role.

5.5.3. Cassette Settings


Access the list of stored Cassette Settings. The list is according to the scan time, with the most recent scan on top. Scroll up and down across the touchscreen to move through the list and press on a row to explore the Cassette Settings details.

The following buttons may be available in this screen:

- Restore:** Restores the factory Cassette Settings list; in case it is an empty list all Cassette Settings will be removed. If this button is available, it appears on the upper right corner of the screen. (See 'Network' section).
- Add:** Displays the screen to load a new Cassette Settings through the barcode reader.
- Select as current (one per item):** Press this button, to set the corresponding Cassette Settings as the current Cassette Settings. In doing so, the selected Cassette Settings will be ready for use and its name will appear on the top bar.

NOTE: The maximum number of Cassette Settings that LFR stores could be set between 1 and 20 (See 'Cassette' section).

To access the details screen of a specific Cassette Settings, press on its corresponding row on the list. When a Cassette Setting is selected, several tabs are available:

- Lines:** The is tab shows the list of test lines included on this Cassette Settings, grouped by strips.
- Batch:** This tab shows the information related to the batch of this Cassette Settings.
- Info:** This tab shows any informative text included in Cassette Settings.

LINES BATCH INFO

Batch ID: 03A9

Expiration Date: 2027-06-30 02:00

Export Date: 2025-08-27 10:46

DELETE

LINES BATCH INFO

Description: Neutral settings to use in IFU documents. Qualitative

DELETE

In addition, the screen also displays a ‘Delete’ button to delete the selected Cassette Settings. When pressed, a confirmation screen appears. Press ‘Confirm’ to accept the deletion of the selected Cassette Setting or ‘Cancel’ or anywhere outside the confirmation screen to cancel the deletion.

Delete

Confirm to delete the current cassette settings?

ⓘ CANCEL

✔ CONFIRM

5.5.4. Audit Trail

1. Operations

This table shows the records of the various operations performed on the LFR. It is possible to reorder the information by clicking on the column headers. Each column can be sorted in ascending or descending order according to the value of its cells.

ADD FILTER

Date/Time ↓	Action	Destination	Type	Source IP	Request	Sample ID	User	UUID
03/17/2026, 12:26	SCAN	DB	PROCESS	lfr	POST	test125455	admin	d47ed613-3c06-4199-9dd9-87ad51f5835e
03/17/2026, 12:23	SCAN	DB	PROCESS	lfr	POST	test54242	admin	cce32bfa-f3ed-4087-b6df-2e94624e2d92
03/17/2026, 12:22	SCAN	DB	PROCESS	lfr	POST	test44555	admin	adad2b65-3118-417f-ab9f-12293c1b45ec
03/17/2026, 12:20	SCAN	DB	PROCESS	lfr	POST	test45155	admin	82418524-8c68-4ea9-acdd-3d4478a51825
03/17/2026, 12:19	SCAN	DB	PROCESS	lfr	POST	test415425	admin	1a596469-c435-4a24-866c-47a63e84f6bd
03/17/2026, 12:18	SCAN	DB	PROCESS	lfr	POST	test42542542	admin	db6148a7-71d4-4077-850c-5d909cffce32
03/17/2026, 12:18	SCAN	DB	PROCESS	lfr	POST	test47242	admin	4b6ee4ce-9b8e-4ff5-8b81-f33b1178c862
03/17/2026, 12:17	SCAN	DB	PROCESS	lfr	POST	test187646	admin	9363d80d-2b23-485f-9e36-8fdb999d4d5c
03/17/2026, 12:14	SCAN	DB	PROCESS	lfr	POST	test12238	admin	d991fbd4-c7c1-4b91-b539-44a4a4c49a9d
03/12/2026, 16:00	SCAN	DB	PROCESS	lfr	POST	dhrdf	admin	53cf5269-fb50-4edb-8146-d613064cd18e

Rows per page: 10

1-10 of 23

1

2

3

NEXT >

To facilitate the search for specific information, the display has a filtering function. To activate it, click on the icon in the upper right corner. Upon doing so, a context menu appears to select the parameter (column) on which you want to apply the filter. Once the parameter is selected, a text field appears in the upper left corner of the screen. Click on this field to open the virtual keyboard, enter the desired value and confirm the display of the filtered list.

2. Audits

This section is primarily used for auditing the registered exportation. It functions as a log file to track remote audits.

Date/Time ↓	Action	Destination	Type	Source IP	Request	Sample ID	User	UUID
28/08/2025 12:00	AUDIT	M2MAPI	HL7	192.168.138.13	GET	N.A.	Authorized Machine	N.A.

1-1 of 1 < >

5.6. Tools Menu

Press the gear symbol on the centre-right of the Top Bar to access the Tools Menu with the following options:

 System Menu

 User Management

 Logout

5.6.1. System Menu

Explained in the ‘System’ section.

5.6.2. User Management

This screen allows password change and is different depending on the current user role – administrator, manager or operator:

1. **As an Administrator or Manager:** A username can be selected from the drop-down list. A new password must then be typed in both the ‘New Password’ and ‘Confirm New Password’ fields to update the password of the selected user.

CHANGE PASSWORD USERS


Username

New Password

Confirm New Password

SAVE

2. **As an Operator:** An operator can change their own password by typing the current password in the 'Current Password' field and typing the new password in both the 'New Password' and the 'Confirm New Password' fields.

CHANGE PASSWORD

Current Password

New Password

Confirm New Password


SAVE

After editing all fields, press the 'Save' button:

- If the LFR displays the green message: '**Password Changed!**', the password is replaced successfully.
- If the LFR displays the red warning: '**Password Not Match!**', the 'New Password' field and the 'Confirm New Password' field are not equal.
- If the LFR displays the red warning: '**Authentication failed!**' the 'Current Password' field is not correct.

5.6.3. User Management – Users (Administrator or Manager only)

Press the 'Users' tab to create, edit, and delete user accounts. This screen is only accessible to users with Administrator or Manager role.

CHANGE PASSWORD USERS

Username*

--

Create new user

aaa

abc

To see available the users press on the field 'Usernames' and a drop-down list is displayed.

1. Create New User

To create a new user, click on 'Create new user' under the 'Users' drop-down list, and fill in the fields available. Only the 'Username', 'Password', and 'User Role' fields are mandatory. Press the 'Save' button and then the 'Confirm' button from the pop-up dialog to store the new user account.

CHANGE PASSWORD **USERS**

Username

Password

Forename

Surname

E-mail

Institution

Department

User Role *

CANCEL **SAVE**

2. Edit or Delete a user

To edit or delete a user account, select the user on the 'Username' drop-down list.

Press each field to edit its value using a virtual keyboard. To save changes, press the 'Save' button and then the 'Confirm' button on the pop-up dialog.

To delete the selected user, press the 'Delete' button and then the 'Confirm' button on the pop-up dialog.

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Usuario *

user

Nombre

user

Apellido

user

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user@user.com

Institución

user

Departamento

user

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
Operator

BORRAR **GUARDAR**

5.6.4. Logout

Click on the Logout button to log out of the current user's session. After logging out, you will be taken to the Login Screen. See the 'Login Screen' section for an explanation of the login process.


5.7. Software update

CAUTION 	Run this procedure under Technical Support advice and with a full battery charge and the LFR plugged to the mains supply.
CAUTION 	Software update packages are only compatible within the same major version. Cross-version updates (e.g., installing 4.x on a 3.x device) are not supported.

To update the Software, follow these steps under technical support supervision:

1. Format a slim USB memory stick in the 'System' tab from the System Menu of the LFR (see 'System' section).
2. Under technical support supervision, download the proper update package from the website into a computer and unzip it.
3. Copy the update package (filename.upd) from the computer to the root of the slim USB memory stick.
4. Switch off the LFR.
5. Plug the memory stick into the USB port of the LFR.
6. Switch on the LFR.
7. A software update confirmation screen will appear. Click 'Accept' to continue:

8. An update screen is displayed. This update screen depends on the Wipe Clean Version of the LFR.
9. Once the update is finished, the screen will display 'Update Success'. Follow the instructions indicated on the screen.

5.8. Browser Generic Operations

Here are some tips to browse the iPeak® reader 4,3" .

- Remove an error message from the screen, by touching the screen outside the error.
- Depending on the selected field, a virtual keyboard is displayed:

For Username(user), Password(user), Sample ID, Comments, Name, Surname, and Department:

a / A	d / D	g / G	j / J	m / M	p / P	s / S	v / V	y / Y	1	4	7
b / B	e / E	h / H	k / K	n / N	q / Q	t / T	w / W	z / Z	2	5	8
c / C	f / F	i / I	l / L	o / O	r / R	u / U	x / X	0	3	6	9

For SSID, Password (WLAN), Email, Institution, Server, HL7 path, PDF path, Server (SFTP and TCP/IP), Username (SFTP) and Password (SFTP):

0 ... 9			a ... z			A ... Z				
space	!	"	#	\$	%	&	'	(	)	=
*	+	,	-	.	/	:	;	<	>	^
?	@	[	]	\	_	'	'	{	}	~

For IP, Mask, and Gateway:

1	4	7	0
2	5	8	.
3	6	9	

For Cassette Settings Version and Port:

1	4	7	0
2	5	8	
3	6	9	

6. TROUBLESHOOTING

6.1. Error and troubleshooting list

Symptom	Probable cause	Recommended action
The unit does not switch ON.	Battery is discharged	Plug the Power supply unit to recharge
The LFR does not ON or OFF	ON/OFF button is damaged	Contact the Technical Service for repair
The unit switches off spontaneously (a message has been displayed for a while before switching off)	Too much time of inactivity has elapsed	Increase the inactivity time in the System Options menu
	Mains cord unplugged and battery level is too low	Plug the mains supply and let the battery fully charge
Red LED is blinking while battery is charging	Incorrect charger is used	Ensure that the charger used is the correct one
	The charger or power cord is not properly plugged	Ensure that all charger wiring is well connected
	The charger or power cord is damaged	Replace it with a new one
	The battery or other electronic component is damaged	Contact the Technical Service for repair
The unit does not perform as expected	Old version of SW	Upgrade the SW
The Wi-Fi SSID of the LFR does not appear on the computer	The wireless signal is out of range	Place the unit closer to the wireless connection
	The Wi-Fi device is damaged	Contact the Technical Service for repair
The Ethernet is not detected	The Ethernet is not connected properly	Ensure that the Ethernet connector is well placed in the socket.
	Damage to Ethernet socket or its internal connection port	Contact the Technical Service for repair
After reading barcode for Cassette Settings or for confirmation, an error message is shown	The barcode type is different from the expected one	If you are loading the Cassette Settings, show the big Data Matrix (configuration barcode). If you are confirming the cassette, show the Small-Data Matrix
The barcode reader light is on, but the barcode cannot be read	The barcode is incorrectly placed in front of the barcode reader	Slightly tilt and move the barcode forward-backward, until the barcode is read

Symptom	Probable cause	Recommended action
	The barcode is damaged or is of the wrong type	Replace the barcode with a new or correct one
	Dirt on the barcode reader window	Clean the barcode reader with an antistatic cloth or a brush
The LFR cannot access the USB memory stick	The USB stick is too wide, hence it is not properly inserted into the socket	Use a slim USB memory stick
	The USB is not properly formatted	Format the USB memory stick using the option in the System options
	The USB is not properly inserted	Correctly insert the USB memory stick
	Damage to USB socket or its internal connection port	Contact the Technical Service for repair
A scan cannot be properly performed and the following message is displayed "Please fit drawer and cassette properly"	Some obstacle is preventing the drawer from being inserted correctly and/or completely	Remove the obstacle
	No drawer has been inserted in the LFR	Insert the cassette in the drawer correctly and completely
	The drawer is damaged	Replace the drawer
When launching a scan, the device shows a message: Normalization error: fiducial not found	The white dots of the drawer are dirty or damaged	Clean or replace the drawer
When displaying results, the Control line says: Absent or Invalid, and no results for test lines are shown	The control line is absent or too weak	Correctly inoculate the cassette. Wait for the prescribed time before launching the scan or use the 'Timed Scan' feature
The screen displays the message "RTC CRITICAL FAILURE Please contact Technical Service"	Real Time Clock battery is damaged	Contact the Technical Service for repair
After performing a scan, the screen displays the error message "Error trying to export results to SFTP server..."	No internet connection	Review network connection and read the entire error message for more options
	The SFTP server is not properly configured	Review the SFTP configuration in the 'Network' tab from System Menu
After performing a scan, the screen displays the error message	No internet connection	Review network connection and read the whole error message for more options

Symptom	Probable cause	Recommended action
"Error trying to export results to TCP/IP server..."	The TCP/IP server is not properly configured	Review the TCP/IP configuration in the 'Network' tab from System Menu
The LFR cannot reset to default settings	The internal button is not properly pressed	Ensure the activation of the button during start-up
	Damage to the Software Wipe internal button	Contact the Technical Service for repair
The device restarts itself	Electrostatic discharge on the touchscreen	Wait until the reboot process ends
The device does not respond	Electrostatic discharge on the touchscreen	Press and hold the power button for 10 seconds. Press the power button again to switch on the device

6.2. Remote Support

To perform remote support following these instructions:

- Under 'System Menu' then Preferences, disable the option 'Inactivity shutdown enable'
NOTE: This is important to ensure the reader does not switch off during the remote session.
- Connect the reader to the internet via Ethernet cable (not Wi-Fi).
- Under 'System Menu' then 'System', press the 'Connect' button of the 'Remote Support' function.
- After clicking 'Connect', a series of numbers will appear on the screen. Send the Technical Support the following:
 - The full number sequence (this corresponds to the port number)
 - A photo of the screen showing these numbers.
- Once connected, do NOT switch off the reader or disconnect it from the network, as this will interrupt the remote access.
- As soon as connection details are received, the technical team will access the system and proceed with the analysis.

7. MAINTENANCE

7.1. Cleaning Procedure

Clean the device and drawer surfaces using a soft, lint-free cloth in accordance with applicable laboratory cleaning practices.

Do not use abrasive materials, excessive moisture, or cleaning methods that may damage the device surfaces or optical components.

CAUTION 	Damage to the device Do not spray inside the LFR or close to the drawer slot.
---	---

CAUTION 	Damage to the drawer Clean the 4 white dots very carefully to avoid damaging them.
---	--

7.2. Service Life and Preventive Maintenance

The iPeak® reader 4,3" does not have a predefined expiration date. However, to ensure continued safe operation and optimal performance over time, the manufacturer establishes a mandatory preventive maintenance program.

The iPeak® reader 4,3" is guaranteed to operate safely for five (5) years from the manufacturing date without requiring maintenance during that period. After this initial period, preventive maintenance must be performed every two (2) years, provided that critical components remain available and unaffected by obsolescence. Maintenance may include, where applicable, functional checks, software updates, cleaning, readjustment, or the replacement of wear-sensitive parts. Use of the iPeak® reader 4,3" beyond the five-year period without performing the required maintenance compromises its conformity and is entirely at the user's own risk. The manufacturer is not liable for performance, safety, or regulatory issues resulting from missed maintenance.

The iPeak® reader 4,3" is intended to be used as part of a broader diagnostic system. It must be integrated and validated by the system developer or test manufacturer. The iPeak® reader 4,3" manufacturer is not responsible for the regulatory compliance or performance of the overall system.

NOTE: This service life and maintenance guidance applies only to the main instrument. It does not affect the separate shelf life, storage conditions, or expiration dates of any accessory, consumables, reagents, control cassettes, or third-party components provided by IUL or its partners, which remain subject to their respective specifications and regulatory markings and do not affect the declared service life of the device.

APPENDIX A ORDERING CODES

Spare parts

	Product	Contents	Cat. No.
	Customized drawer	1 unit	9000xxxxx
	Colorimetric self-adjustment light drawer (Standard version)	1 unit	900013732
	External power supply (3-poles)	1 unit	900020088
	External power supply (2-poles)	1 unit	900020089
	2-poles power cord (EU)	1 unit	900020090
	2-poles power cord (US / JP)	1 unit	900020080
	2-poles power cord (UK)	1 unit	90002264
	3-poles power cord (EU)	1 unit	900020077

	Product	Contents	Cat. No.
	3-poles power cord (UK)	1 unit	900020079

Optional Accessories

Product	Contents	Cat. No.
Ticket DOT Matrix printer (serial interface, 40 columns)	1 unit	90004875
Ticket Thermal Printer (serial interface)	1 unit	90005048
Colorimetric UQC drawer (Standard version)	1 unit	900012106

APPENDIX B TECHNICAL DATA

The manufacturer reserves the right to change specifications at any time.

Specifications

Process time	from 15 seconds
Camera sensor	8 Mpixel
Touch screen	4,3"
Operating system	Linux
Coefficient of variation (CV)	< 5% ¹

Connectivity

Ports	Wi-Fi, Ethernet, and slim USB
Export results	HL7 v2.8 compliant. Connections https with SSL TLSv1 TLSv1.1 TLSv1.2 TLSv1.3. Authentication: transport level 2
Audit trail	HL7 v2.8 compliant. Connections https with SSL TLSv1 TLSv1.1 TLSv1.2 TLSv1.3. Authentication: transport level 2
SFTP	LFR can export to an external SFTP server
TCP/IP	LFR can export to an external TCP/IP server

Operating Conditions

Power range	12V +/-5% (it is supplied with an external power supply unit that can be used in 100V-240V)
Frequency range	D.C. (it is supplied with an external power supply unit that can be used in 50Hz-60Hz)
Maximum power	40 W

¹ In lines with an intensity around the center of the dynamic range.

Overvoltage category	II
Ambient temperature	20 - 25 °C / 68 – 77 °F
Relative humidity	10 - 75 % (non-condensing)
Altitude	Up to 2000 meters (6500ft.)
Place of operation	For indoor use only
Pollution level	2

Transportation conditions

Temperature	5 – 40 °C / 41 – 104 °F In manufacturer's packaging
Relative humidity	Maximum 75 % (non-condensing)

Storage conditions

Temperature	5 - 40°C / 41 - 104°F In manufacturer's packaging
Relative humidity	Maximum 75 % (non-condensing)

Dimensions and Weight

Dimensions (LxWxH)	12 x 15 x 14 cm / 4.6 x 5.8 x 5.5 in
Weight	0.6 Kg / 1.3 lb

APPENDIX C WARRANTY

The iPeak® reader 4,3" has a 12-month warranty period.

The warranty is void when misuse of the equipment can be proved. Damage or faults caused by impacts, chemical or corrosive products, liquids, damp, or other external factors, such as radiation, fire, or inadequate transport are not included.

In addition, the warranty will not apply if the equipment has been handled, repaired, or modified by not qualified or specifically designated personnel.

As part of this warranty, all shipping costs to the selected technical service are the responsibility of the Buyer.

APPENDIX D WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE)

This appendix provides information about the disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

APPENDIX E EU DECLARATION OF CONFORMITY

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Single Registration Number (SRN): ES-MF-000000861

Declares that the Product Name:

iPeak® reader 4,3"

Catalog Number:

REF: 900033000 ¹

Basic UDI-DI:

8437022431LFR8L

Intended Purpose:

iPeak® reader 4,3" is an instrument specifically to be used to provide quantitative, semi-quantitative, and/or qualitative in-vitro determination of photometric immunochromatographic test defined and marketed by IUL's partners.

Classification:

According to Rule 5(b) of Annex VIII, the product is classified as CLASS A.

Statements:

- EU declaration of conformity is issued under the sole responsibility of the manufacturer.
- The product that is covered by the present declaration is in conformity with the Regulation (EU) [2017/746](#) of the European Parliament and of the council of 5 April 2017 on in vitro diagnostic medical devices.

Barcelona, November the 28th 2022

A blue ink signature of Vicenç Font, written in a cursive style.

Vicenç Font
(IVD Technical Manager)

A blue ink signature of David Guerrero, written in a cursive style.

David Guerrero
(Quality Manager)

¹ Generic reference, current version of this EU Declaration of conformity available in IUL internal document 50008290.

APPENDIX F FCC STATEMENT

FCC Declaration of Conformity

Responsible party:

IUL, S.A.

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08030 Barcelona

España

This device:

- Contains TX FCC ID:2ABCB-RPI4B
- Contains IC: 20953-RPI4B
- Complies with Part 15 of FCC Rules, Operation is Subject to the following two conditions:
 - (1) This device may not cause harmful interference, and
 - (2) This device must accept any interference received including interference that causes undesired operation.

Any changes or modifications to the equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the device.

APPENDIX G UK DECLARATION OF CONFORMITY

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UK Responsible Person:

CMC Medical Devices Ltd. – Office 32 19-21 Crawford Street, W1H 1PJ London, United Kingdom

Declares the following product:

iPeak® reader 4,3" (REF: 900033000¹)

Intended Purpose:

The iPeak® reader 4,3" product is an instrument specifically to be used to provide quantitative, semi-quantitative, and/or qualitative in-vitro determination of photometric immunochromatographic test defined and marketed by IUL's partners

Classification:

NOT LIST A or B, ARTICLE 9 – GENERAL IVD DEVICES

Conformity assessment route:

IN VITRO MEDICAL DEVICE DIRECTIVE 98/79/EC, ANNEX III (SELF DECLARATION)

Statements:

- This declaration of conformity is issued under the sole responsibility of the manufacturer
- The product that is covered by the present declaration is in conformity with the UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)

UK designated standards:


- EN ISO 13485:2016 & EN ISO 13485:2016/AC:2018 - Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 61010-2-101:2002 - Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment
- EN 61326-2-6:2006 - Electrical equipment for measurement, control, and laboratory use - EMC requirements - Part 2- 6: Particular requirements - In vitro diagnostic (IVD) medical equipment
- EN 62304:2006 - Medical device software - Software life-cycle processes.

Barcelona, February the 14th 2024

**UK
CA**

Vicenç Font
(IVD Technical Manager)

David Guerrero
(Quality Manager)

¹ Generic reference, current version of this UK Declaration of conformity available in IUL internal document 50008514.

APPENDIX H ROHS STATEMENT

The following information has been made available to comply with The Restriction of Hazardous Substances Directive, (RoHS2 & RoHS3), short for Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

APPENDIX I APPLIED STANDARDS

The iPeak® reader 4,3" is in conformity with the following standards:

- ISO 13485:2016 Medical devices - Quality management systems - Requirements for regulatory purposes
- CFR21 PART 820 Quality System Regulation
- EN 61010-1:2010 Safety requirements for electrical equipment for measurement, control, and laboratory use -- Part 1: General requirements
- EN 61010-2-101:2017 Safety requirements for electrical equipment for measurement, control, and laboratory use — Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment
- EN 60601-1:2006+AC:2010+A1:2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- EN 61326-1:2013 Electrical equipment for measurement, control, and laboratory use - EMC requirements - Part 1: General requirements
- EN 61326-2-6:2013 Electrical equipment for measurement, control, and laboratory use — EMC requirements — Part 2-6: Particular requirements — In vitro diagnostic (IVD) medical equipment (IEC 61326-2-6:2012)
- IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
- ISO 14971:2019 Medical devices – Application of risk management to medical devices
- IEC 62304:2006 & IEC 62304:2006/A1:2015 Medical device software – Software life-cycle processes
- ISO 15223-1:2021 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements.
- CAN/CSA C22.2 No. 61010-1: Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements - Third Edition; Update No. 1: July 2015; Update No. 2: April 2016; Amd 1: November 2018
- CAN/CSA C22.2 No. 61010-2-101-19 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-101 Particular requirements for in vitro diagnostic (IVD) medical equipment.
- UL 61010-1 UL Standard for Safety Electrical Equipment For Measurement, Control, and Laboratory Use; Part 1: General Requirements - Third Edition; Including Revisions through November 21, 2018
- UL 61010-2-101-19 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-101 Particular requirements for in vitro diagnostic (IVD) medical equipment

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
The login screen features a blue circular icon with a white padlock at the top. Below it, the word 'Bienvenido' is displayed in a light blue font. There are two input fields: the first is labeled 'Usuario' in red, followed by a red line and the word 'Requerido' in red; the second is labeled 'Contraseña' in red, followed by a red line and the word 'Requerido' in red. At the bottom, there is a grey button with the text 'ACEPTAR' in blue capital letters.

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IUL Inc.
IUL GmbH